

Dental Unit (Economic Model)

MST-EM215A



Standard Accessories:

Two-way or Four-way high speed air turbine H.p:	2 sets
Two-way or Four-way low speed air motor H.p:	1 set
Three way spring:	1 set
Water suction and saliva ejector:	1 set
Operation light:	1 set
X-film viewer:	1 set
Water purified system:	1 set
Dental stool:	1 piece

Specifications:

Supply Pressure: ~220V, 50Hz
 Air pressure: >0.55MPa
 Water Pressure: 0.2MPa~0.4MPa

Dental Unit (Economic Model)

MST-EM216



Standard Accessories:

4 hole high speed air turbine hand piece:	2 sets
4 hole low speed air turbine hand piece:	1 set
3-way syringe(cold/hot water):	2 sets
Powerful suction device:	1 set
Salivea ejector:	1 set
Dentists seat:	1 set
Dental light:	1 set
Warm water system:	1 set
Water supply system:	1 set
Rotatable ceramic spittoon:	1 set

Specifications:

Supply Pressure: ~220V,50Hz
 Air pressure: >0.55MPa
 Water Pressure: 0.2MPa~0.4MPa



Dental Unit (Supreme Model)

MST-SM9002

Factory Equipment	
Dc driven chair:	1 set
Motor dental chair control system:	1 set
With three memory	
X-ray viewer:	1 set
Auto spittoon flush and cup filler control system:	1 set
Clean water system:	1 set
Three way syringe:	2 sets
Saliva ejector:	1 set
High suction:	1 set
Dental light:	1 set
Dental stool:	1 set

Options

- High-speed hand piece(4-holes)
- Low-speed hand piece(4-holes)
- Intra-oral camera
- Ultrasonic scaler
- Light curing
- Dental air compressor

Specifications

- Power voltage:220v 50Hz
- Motor voltage 24v
- Water pressure 0.2-0.4mpa
- Air pressure:0.5-0.8mpa

Dental Unit (Supreme Model)

MST-SM9004



Factory Equipment	
Dc driven chair:	1 set
Motor dental chair control system:	1 set
With three memory	
X-ray viewer:	1 set
Auto spittoon flush and cup filler control system:	1 set
Clean water system:	1 set
Three way syringe:	2 sets
Saliva ejector:	1 set
High suction:	1 set
Dental light:	1 set
Dental stool:	1 set
Multifunction kick switch:	1 set

Options

- High-speed hand piece(4-holes)
- Low-speed hand piece(4-holes)
- Intra-oral camera
- Ultrasonic scaler
- Light curing
- Dental air compressor

Specifications

- Power voltage 220v 50Hz
- Motor voltage:24v
- Water pressure:0.2-0.4mpa
- Air pressure 0.5-0.8mpa

Dental Unit (luxurious model)

MST-LM398HA

Standard Equipment:	
AL-318 endoscope system:	1 set
4 hole high speed air turbine handpiece:	2 set
4 hole low speed air turbine handpiece:	1 set
Curing light (built-in):	1 set
Ultrasonic scaler (built-in):	1 set
Tempered glass spittoon:	1 set
Auto spittoon flush and cup filler control system:	1 set
Motor dental chair with 3 adjustable positions:	1 set
Danish 24V noiseless DC motor chair:	1 set
Extendable cushion and backrest:	1 set
High-grade articulated headrest:	1 set
Operating light:	1 set
Saliva ejector with high and low suction:	1 set
Instrument tray with air brake:	1 set
3-way syringe (hot/cold water):	2 set
Full-functional foot controller:	1 set
Water tank (air pressurized):	1 set
Operating stool 338C:	1 set

Specifications:
 Power voltage 220V / 50Hz
 Motor voltage 24V
 Water pressure 0.2MPa-0.4MPa
 Air pressure 0.5MPa-0.8MPa



Dental Unit (luxurious model)

MST-LM398Sanor'e

Standard Equipment:	
4 hole high speed air turbine handpiece:	2 set
4 hole low speed air turbine handpiece:	1 set
Imported electromagnetism valve:	1 set
Rotatable tempered glass spittoon (hitting obstruction safety-equipment):	1 set
Auto spittoon flush and cup filler control system:	1 set
Motor dental chair with 3 adjustable positions:	1 set
Danish 24V noiseless DC motor chair:	1 set
Foldaway patient chair synchronized movement of seat and backrest:	1 set
Connected ground type cabinet:	1 set
High-grade articulated headrest:	1 set
Operating light:	1 set
Saliva ejector with high and low suction:	1 set
Noise elimination and filtration system of high power suction:	1 set
Instrument tray with air brake:	1 set
3-way syringe (hot/cold water):	2 set
Full-functional foot controller:	1 set
Imported 24V X-film viewer:	1 set
Water tank (air pressurized):	1 set
Operating stool 338C:	1 set

Options:

- Scaler (built-in)
- LED curing light (built-in)
- Endoscope system

Specifications:

- Power voltage: 220V / 50Hz
- Motor voltage: 24V
- Water pressure: 0.2MPa-0.4MPa
- Air pressure: 0.5MPa-0.8MPa