



HX-D01 手摇升降急诊床车	HX-D01 Manual lift emergency bed
· 规格: (mm) 2010×570×550/780 床面: 1900×570	SPECIFICATIONS: (mm) 2010×570×550/780 Bed board:1900X570
· 功能: 单摇	FUNCTION: single shakes
· 材质: 不锈钢	MATERIAL: stainless steel
· 标准配置: 输液架一支; 海绵垫一块	STANDARD CONFIGURATION: a transfusion pole, a sponge cushion



HX-D04 不锈钢对接过床车	HX-D04 Stainless steel connecting stretcher
· 规格: (mm) 3600×570×550/780 床面: 1900×570	SPECIFICATIONS: (mm) 3600×570×550/780 Bed board:1900X570
· 材质: 输液架一支	MATERIAL: a transfusion pole,



HX-D02 不锈钢两大两小轮担架车 (病人运输车)	HX-D02 Stainless steel stretcher trolley with two big wheels and two small wheels
· 规格: (mm) 2000×700×800 床面: 1900×600	SPECIFICATIONS: (mm) 2000×700×800 Bed board:1900X600
· 材质: 不锈钢	MATERIAL: Stainless steel



HX-D05 横式过床车	HX-D05 Horizontal stretcher
· 规格: (mm) 2040×660×550/800 床面: 1900×570	SPECIFICATIONS: (mm) 2040×660×550/800 Bed board:1900X570
· 材质: 不锈钢; 配输液杆一支	MATERIAL: Stainless steel; a transfusion pole



HX-D03 四小轮担架车 (折式出口装)	HX-D03 Stretcher trolley with four small wheels
· 规格: (mm) 2000×680×800 床面: 1900×600	SPECIFICATIONS: (mm) 2000×680×800 Bed board:1900X600
· 材质: 不锈钢	MATERIAL: Stainless steel
· 标准配置: 输液架一支; 5寸万向轮 (对角刹车)	STANDARD CONFIGURATION: a transfusion pole, a sponge cushion



HX-D06 塑钢手摇升降急诊床车	HX-D06 Steel jet molding manual lift emergency bed
· 规格: (mm) 1900×610×550/860	SPECIFICATIONS: (mm) 1900×610×550/860
· 功能: 背部倾斜80°±5°; 整体升降310mm.	FUNCTION: back rest 80°±5°; the overall take-off and landing 310mm
· 材质: 床架碳钢材料制作, 采用抗菌环保粉体涂装, 床面为ABS面板;	MATERIAL: carbon steel tube bed frame, most up-to-date double anti-bacterial and anti-rust powder coating treatment; ABS plastic cover base.
· 标准配置: 6寸豪华中控轮; ABS欧式护栏两片; 输液杆一支。	STANDARD CONFIGURATION: 6" luxurt castor ;two ABS European guard rails a transfusion pole